

	GoTalk Pocket	GoTalk 4+	GoTalk 9+	GoTalk 20+	Super Talker	GoTalk Express	Tobii S32
Größe/ Gewicht	14,5 x 9,5 cm <200 g	23 x 31 x 30 cm 340 g	23 x 30 x 3 cm 650 g	20 x 30 x 3 cm 650 g	305 x 229 x 51 mm 900 g	32 x 25,5 x 4,5 cm 1,1 kg	31 x 25 x 3 cm 862 g
Felderanzahl	6	4 + 2	9 + 3	20 + 5	1, 2, 4 oder 8	32	1, 2, 4, 8, 16 oder 32
Sprach- speicher	5 min	5,5 min	9 min	15 min	16 min	26,6 min	> 60 h
Ebenen	5	4	5	5	8	5	192
Sonstige Funktionen	-	-	-	-	Tasteranschlüsse für jedes Feld  2 Ausgangsbuchsen	Scanning  Expressfunktion	IR-Schnittstelle für Umfeldsteuerung Scanning optional Sequenzen visuelle Szenen
Software	Nicht im Lieferumfang	Nicht im Lieferumfang	Nicht im Lieferumfang	Nicht im Lieferumfang	Nicht im Lieferumfang	Nicht im Lieferumfang	Symbol Mate enthalten
							